

Objektivadapter für Mikroskope M3 / MS5 / MZ6 / MZ7.5

Objective Adapters for Microscopes M3 / MS5 / MZ6 / MZ7.5

solino™




Diese Adapter können mit Mikroskopen M50x1 Objektivgewinde benutzt werden. Z.B.:

- Wild Leica M3B, M3C, M3Z
- MS5, MZ6, MZ7,5

These adapters can be used with all Microscopes using the M50x1 objective thread. This includes:


- *Wild Leica M3B, M3C, M3Z*
- *MS5, MZ6, MZ7,5*

Microscopes M3/ MS5 / MZ6 / MZ7.5
M50x1



MZ-Objectives
M50x1

043-300050	Adapter OD M50x1 / ID M60x1
043-300062	Adapter OD M50x1 / ID M65x1.5
043-300077	Extension 21.5 mm for M50x1, OD 58 mm



	A	B	C	D
043-300050	M60x1	M50x1	12.5 mm	62 mm
043-300062	M65x1.5	M50x1	12.5 mm	69 mm
043-300077	M50x1	M50x1	21.5 mm	58 mm