

Objektivadapter für Mikroskope M125 / M165 / M205 Objective Adapters for Microscopes M125 / M165 / M205

solino™




Diese Adapter können mit Mikroskopen mit M65x1,5 Objektivgewinde benutzt werden. Z.B.:


- Leica M125 / M165 / M205
- Leica MZ 12.5 / MZ 16
- Leica SOM-Microscopes

These adapters can be used with all Microscopes using the M65x1.5 objective thread. This includes:

- *Leica M125 / M165 / M205*
- *Leica MZ 12.5 / MZ 16*
- *Leica SOM-Microscopes*



043-300024	Adapter OD M65x1.5 / ID M50x1
10 450 296	Adapter OD M65x1.5 / ID M60x1
043-300078	Extension 15 mm M65x1.5, OD 69 mm
043-300076	Extension 30 mm M65x1.5, OD 69 mm



	A	B	C	D
043-300024	M50x1	M65x1.5	6 mm	69 mm
10 450 296	M60x1	M65x1.5		
043-300078	M65x1.5	M65x1.5	15 mm	69 mm
043-300076	M65x1.5	M65x1.5	30 mm	69 mm