

Mobile Floor Stand Advanced, Mount for Extensions

Fahrbares Bodenstativ Advanced m. Auslegeraufnahme




Technical Features

Article number 043-603020

- Tripo stand with lockable wheels
- Height stepless adjustable from 1280 mm to 2130 mm
- Arm 700 mm with optional plate for cold light source Schott series or tools

Technical Drawing



Configuration

043-603008
Height Adjustable Swivel Arm; 5/8"-Hole

043-603008-132
Height Adjustable Swivel Arm; Column D=32 mm

043-603008-125
Height Adjustable Swivel Arm; Column D=25 mm

043-603020
Mobile Floor Stand Advanced, Mount for Extension Arms

043-530114
Universal Lightsource mounting plate,
210x210 mm with adjustable velcro strap

