



Focus Drive C, Coarse

Fokussiertrieb C, Grob



Technical Features

Article number 043-610004	
Focus Drive C	coarse
Travel range	57 mm
Travel per knob rotation	21 mm
Load	max. 10 kg
Accessories for attachment to stands and walls	
043-610904-001	Adapter plate to wall and Track Stand
043-610904-320-068	Clamp to Ø 32 mm column and Post Stand
Accessories for carrying microscopes	
050-502076126W12	LED-ring-light incl. clamp with Ø 76mm
043-610920	Microscope carrier with Ø 76 mm
043-610914	Clamp mini with Ø 35 mm for TV-microscopes



Technical Drawing

