

LEDrive

Focus drive C, coarse w/integ. LED ring and 5/8" pin adapter



High Class Microscope Accessories at Opto


LEDrive with integrated ring light and 5/8" interface is specially geared to the efficient and ergonomic use of a stereo microscope, in order to make the handling even easier and more flexible.

LEDrive is part of the Opto **Microscope Accessories** for stereomicroscopes of leading microscope manufacturers such as Leica, Nikon, Olympus, etc.




Technical Features

Article number 043-643045	
Focus drive C	coarse
Travel	57 mm
Load	max. 9 kg
Universal scope carrier	Ø 76 mm
Ring light, integrated into carrier	LED, color white approx. 6800 K
Illuminance	22500 lx, at working distance 245 mm
	40000 lx, at working distance 160 mm
	61000 lx, at working distance 110 mm
Power consumption	12 V / 10 W
Cable length	2 m
Stand adapter, pin	5/8"
Accessories	
050-70011220TV	Power supply 12 V / 25 W
035-8645BN	Ecoline Stereo Microscope 645BN
043-603003	Mobile Floor Stand



Accessories

Mobile Floor Stand Configuration for flexible surface and component inspection.



Technical Drawing

