


Heavy Duty Focus Drive OP, Coarse/Fine

Heavy-Duty-Fokustrieb OP, Grob/Fein



Technical Features

Article number 100-FB-OP

Focus Drive	coarse/fine
Resolution	< 1 μm
Repeatability	1 μm
Backlash	< 5 μm
Travel range	50 mm
Travel per knob rotation, coarse	4 mm
Travel per knob rotation, fine	0.1 mm
Load	max. 50 kg with pre-adjustment for 20 kg
Weight	2 kg

Accessories for attachment to stands and walls

045-600157-01	Adapter to Column for Focus Drive OP
---------------	--------------------------------------



Technical Drawing

