


## Post Stand - Compact and Modular Microscope Stand

### Kompakte Stativlösung mit 32 mm Rundsäule






## Technical Features

Post Stand – Compact and Modular Microscope Stand for precise measurement and microscopic analysis

Clamp diameter	Ø 35 mm
Vertical travel of Focus Drive C	57 mm
Stand height	510 mm
Column diameter	Ø 32 mm
Travel range of manual stage	75 x 55 mm
Alu plate dimension	300 x 300 x 12 mm
Weight	9 kg

### Accessories

100-JC2-66	Jointed Coupler, 2 links, M6-M6, L=132 mm
------------	---



### Technical Accessories

