


Tilting Table

Ergonomic work for
360° inspection tasks

Workplace with stable
tripod and microscope



Tilting Table - Vacuum with magnetic fixings and hand rests



Perfect Flexibility

Three dimensional objects often benefit from a specific angle of view in order to reveal key details. With our Tilting Table, it is possible to position your microscope securely at an angle or to tilt and fix your samples at the desired positions.

The entire system from granite slab, stable tripod and microscope mounting adds up to a tilt-proof and ergonomic microscope workstation. Its objects and the microscope can be aligned and fixed to each other at a certain angle.

Key Features

- Rock solid stand for many types of inspection
- Inclined viewing
- Easy positioning of objects


Your Benefit

- Highly ergonomic inspection stand
- Quick and easy full 360° inspection
- Compatible with all major microscope brands
- Perfect for Electronic Board Repair



Tilting table | specification

Article no.:	043-631006	043-631003	043-631005
Diameter	ø 300 mm	ø 450 mm	ø 400 mm
Plate	steel	steel	granite
Component fixation	magnetic fixing	magnetic fixing	non-slip rubber ESD surface
Fixing the tilting position	Vacuum adjustment with foot swich	Vacuum adjustment with foot swich	own weight
Weight	ca. 5 kg	7.8 kg (including support plates)	ca. 9 kg
Accessories	2 Magnetic hand support plates 8 Conductor board holding magnets	2 Magnetic hand support plates 8 Conductor board holding magnets	



Version 043-631005
with granite slab and
non-slip rubber ESD surface