




Stage micrometer, 50mm scale with 0.1mm and 0.01mm Objektmikrometer, 50mm in 0,1 und 0,01 mm



Technical Features

Article number 010-310345

- glass measuring scale 50 mm in 0.1 mm and 0.01 mm
- perfect for calibration purposes

Specification:	Scale length	Tolerance per Intervall	Tolerance entire length
	to 5 mm	± 0.001 mm	± 0.003 mm
	from 5 to 20 mm	± 0.002 mm	± 0.005 mm
	from 20 to 50 mm	± 0.003 mm	± 0.010 mm
	Reference temperature: 20° C		

010-310345-2 Stage micrometer incl. DKD Certification