

# PC example

Cleanalyzer

<b>Operating system:</b>	Windows 7 Professional (64 bit)
<b>CPU:</b>	min. Intel®Core i7-4xxx
<b>Main storage:</b>	>=16 GB DDR3-RAM, 1600 MHz
<b>1st Hard disk:</b>	SSD 120 GB
<b>2nd Hard disk:</b>	1 TB, 7200 rpm
<b>Screen resolution:</b>	1920 x 1080 Pixel (Full HD)
<b>Font size:</b>	Small fonts
<b>Graphic card:</b>	<ul style="list-style-type: none"><li>• &gt;=2 GB own RAM and own graphic processor</li><li>• incl. backbuffer</li><li>• VGA/ DVI/ HDMI/ DisplayPort</li></ul>
<b>PCIe</b>	2 free PCIe-Ports
<b>LAN/WLAN</b>	Internet access necessary for support / remote servicing
<b>.NET-Version:</b>	Min .NET Framework 4.0
<b>Operating system rights:</b>	Local administration rights (necessary only for installation)
<b>Positioner interface</b>	USB2.0 or USB3.0
<b>Camera interface:</b>	USB3.0
<b>Report software:</b>	MS Excel

Specifications are subject to change without prior notice.

