


Oblique Viewing Module

Schrägsichtmodul

Module à vue oblique



Technical Description

Article No. 043-112108

- Switchable between oblique and straight view
- Light point in oblique view for finding actual position, battery powered
- Working distance in oblique view 100 mm, in straight view 115 mm
- Integrated 0.3x objective
- Viewing angle 30°
- LED-ring light with power supply 100-240V / 60W
- Extension column-optical axis min. 140 mm
- Fits to microscope series (Leica M50, M60, M80 / Zeiss Discovery series / Olympus SZX series / Nikon SMZ 800/1000)
- Power cord EU included
- Weight: 1.6 kg including adapter

Technical Drawing

