

Long working distance Objective M600



Technical Features

Article number	Working distance	Magnification
043-213001	800 mm	0.13 x
043-213002	500 mm	0.2 x
043-213003	1000 mm	0.1 x
043-213005	1200 mm	0.08 x
Working distance equates approx. to focal length		
Objective thread	M65 x 1.5 mm	
Outer diameter	70 mm	
Height	12.5 mm	
Front faced inner thread	M62 x 1.5 mm	
Ringlight adapter 043-300090	66 mm	
Adapter M50-M60-M80-DMS	10450243	



Technical Drawing

