

DualStation

Stereo- and Inverse Microscope

Transmitted e-mirror light base with Leica M80



Technical Features

Article number, base	045-100181
Transmitted light base with Stereo Microscope	Leica M80, Binocular tube and 10x Oculars *
Focus-objective magnification	1x fine focus travel range
Working distance	90 – 110 mm
Field of view with Leica M80, 10x/23 Ocular	2 – 20 mm
Illumination, Transmitted DF	LED, e-Mirror Technology with adjustable brightness and contrast
Invers microscope with included camera	5 MP
Camera connector	USB 3
Inverse microscope magnification	400x on 17" monitor
Resolution	800 nm
Focus	adjustable
Field of view	0.8 x 0.6 mm
Illumination	LED Köhler-illumination, adjustable brightness and high contrast
Surface material	stainless steel
Ready for RFID-reader integration	

* large wedge not included like in the image

Technical Drawing

