

## Dual-Camera Coupler 1x for Leica Ergo Tube HC

Dual-Kamera-Koppler 1x für Leica Ergo Tubus HC



### Technical Features

Article number 045-200208

Dual-Camera Coupler	optimized for VIS splitting ratio 50/50 others on request
Dimensions	69 x 58 x 40 mm
Magnification	1x
Camera interface	C-mount
Max. sensor size	1"
Adjustment range C-mount	$\pm 2$ mm
Aperture diameter C-mount	18 mm

## Technical Drawing

