

Protection Glass Holder D80 mm with Adapter to LED5000

Schutzglashalter an Objektiv D80 mm mit Aufnahme für LED5000



Technical Features

Article number 045-300205

Aperture diameter	76 mm
Clamp diameter	80 mm
Height	30 mm
Protection glass thickness	1 mm

Adapter to ring light Leica LED 5000 and compatible to stereo objectives with diameter 80 mm.

Accessories

043-320027	Protection glass D = 80 mm / set of 5 single pieces
------------	---

Please contact Opto sales team if you are planning to use it with apochromatic lenses

Technical Drawing

