


Magnetic Coupler – Change Micro Objectives Easily

Magnetkopplung für Mikroobjektive



Technical Features

Magnetic Coupler to change micro objectives easily

Outer diameter	35 mm
Extension length	22 mm

Components top part

045-600158-22	M25x0.75 thread outside
045-600158-21	M26x0.7 thread outside

Components bottom part

045-600158-31	M25x0.75 thread inside
045-600158-32	M26x0.7 thread inside
045-600158-33	RMS thread inside

Technical Drawing

