

# PxF Camera-Support UI-xx Module


050-911104

## Product features

Designed to display live images of IDS-cameras or record and save video images

## Comprehensive functionality

- Setting an optimal live image (+ Brightness, + Exposure time, + Contrast, + RGB, + reset to standard values, Gamma, white balance)
- Increasing the frame rate by limiting the image size (section, binning, subsampling)
- Comprehensive range of pre-processing tools
  - Routine generator (Save, load, run of complex image manipulations)
  - Alteration of brightness and contrast
  - Adjust size and rotation
  - Filter
  - Morphological operations
  - Logical and arithmetic operations
  - Euclidean distance transformations
  - Watershed algorithms



04/2017 - Specifications are subject to change without prior notice.