

PxF Manu-Meter Module

050-911107

Product features

This software module allows an interactive measurement of images. It makes no claims to the content and the quality of the images to be analyzed. The advantage of this module is its universal applicability. All configurations can be stored in profiles and accessed at a later time.

Interactive measurement of

- lengths
- angles
- surfaces
- circles, ellipses
- etc.

Tools for image documentation


- scale
- arrow
- text

Virtual ocular standard overlay grids

- EN 10247_2007
- ASTM E 45_2005
- ISO 4967_1998
- JIS G 0555_2005
- NF A-04-106

Reporting tools

- MS EXCEL
- dhs image database



04/2017 - Specifications are subject to change without prior notice.