

# PxF Stitch-Control Module

050-911118

## Product features

A tool for automatic and interactive merging of lined up images into one image. It can be combined individual image segments in both the horizontal direction (x) and vertical direction (y).

## Functionality:

- As a result, a panoramic image will be generated like a single shot
- Assembling panoramas from any number of individual images, only limited on the memory capacity
- Smooth transition selectable
- Substitution for high-resolution camera
- User interface for free positioning, shifting and overlapping of sub-images
- All application configurations can be stored in profiles



04/2017 - Specifications are subject to change without prior notice.